


Dharmacon™ siRNA selection guide

We offer a wide selection of pre-designed siRNA product lines and formats. This quick selection guide will help to determine the best option for your particular needs.



Product formats

SMARTpool

A mixture of 4 siRNA provided as a single reagent; providing advantages in both potency and specificity.

Set of 4

A convenient option for purchasing aliquots of all 4 individual siRNAs targeting a single gene.

Individual siRNAs

Select 1, 2, 3 or 4 individual siRNAs per gene.

Dharmacon™ ON-TARGETplus™ and siGENOME™ SMARTpool™ reagents and 3 out of 4 siRNAs are guaranteed to silence mRNA by at least 75%. See website for details.

horizondiscovery.com

For more information

To find the contact information in your country for your technology of interest, please visit us at horizondiscovery.com/contact-us

Horizon Discovery, 8100 Cambridge Research Park, Waterbeach, Cambridge, CB25 9TL, United Kingdom